



# Easy 200

Edition: 11.24  
Valid until: 11.25

Deliver to:  
Name:  
Street:  
Town: Postal code:  
E-Mail:  
Tel: Fax:

Deliver to:  
Name:  
Street:  
Town: Postal code:  
E-Mail:  
Tel: Fax:

Order date:  
Purchase order:  
Number required:

Marked for:  
Sunrise order no:

☐ Order ☐ Calculation ☐ Built 4 me

• = Standard ☐ = Cost option  
○ = No cost option

Reimbursement code:

Art. no.: 71100201

Max. user weight:

140 kg

USER WEIGHT

kg

Mandatory information to ensure correct configuration



Base Price

## CLIENT MEASUREMENTS (all measurements without cushion)

SEAT WIDTH (measured in mm from outside frame tube to outside frame tube)

320	340	360	380	400	420	440	460	480	500	520	540
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>

SEAT DEPTH (mm)

340	360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>

No upcharge in combination with SW 520 & 540 mm

FRONT SEAT HEIGHT SELECTION: In dependance of front castor size

Front seat height in mm

Castors	Forks	400	410	420	430	440	450	460	470	480	490	500	510	520	530	540
4"	111x32 // 3 holes (1)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	111x45 // 3 holes (1)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	123x45 // 3 holes (2)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	138x45 // 3 holes (2)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5"	111x32 // 3 holes (1)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	111x45 // 3 holes (1)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	123x45 // 3 holes (2)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	138x45 // 3 holes (2)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
6"	166x45 // 4 holes (2)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	111x45 // 3 holes (1)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	123x45 // 3 holes (2)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	138x45 // 3 holes (2)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
7"	166x45 // 4 holes (2)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

REAR SEAT HEIGHT SELECTION: in dependance of rear wheel size

The rear seat height can not be above the selected front seat height and only be between 10 mm and 110 mm below.

Rear Seat Height (mm)

Rear Wheels	400	410	420	430	440	450	460	470	480	490	500	510	520
22"	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
24"	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
25"	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
26"	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

- (1) Made of Carbotechnology  
(2) Compatible with Power Add-On devices

For more information  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
© Sunrise Medical 2025

**Please follow this measuring method:**

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

**Seat width (SW):**

- a) *For the user:* Allow space between sideguard and hip on each side of the chair (approx. 10 - 20 mm altogether)
- b) *On the wheelchair:* From outside frame tube to outside frame tube.

**Seat depth (SD):**

- a) *For the user:* When sitting upright: leave space for 2 - 3 fingers between front edge of the seat upholstery and back of knee (approx. 30 - 60 mm). Care must be taken with adj. upholstery.
- b) *On the wheelchair:* From back tube to front of seat upholstery.

**Front seat height:**

- a) *For the user:* When sitting with 90° kneeangle: From back of knee to bottom of shoe sole.
- b) *On the wheelchair:* From front of seat tube to the ground . (min. 30 mm longer than user measurement) Remember to include the cushion height.

**Rear seat height:**


- b) *On the wheelchair:* From back of seat tube to the ground.  
**Note:** Max. seat height difference rear to front is 130 mm

**Back height:**

- a) *For the user:* From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) *On the wheelchair:* From back of seat tube to upper edge of backrest upholstery.  
**Note:** The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably

**Back angle:**

**Important:** The back height depends on the degree of disability. 0° is the angle 90° to the ground. Negative angles are <90° and positive angles are >90°.



Setting of the chair above is: SD 400 mm and 24" wheel

**Technical Data Overview:**

**Overall dimensions:**

Max overall height: 800 mm  
 Max overall length: 970 mm

**Overall width with:**

For 0° camber: SW + 200 mm  
 For 2° camber: SW + 220 mm  
 For 4° camber: SW + 280 mm

**24" rear wheel**

**Weights:**

Minimum weight:  
 approx. 13 kg

**Folded Width:**

With rear wheels: approx. 300 mm  
 Without rear wheels: approx. 250 mm

Possible seat height settings. **IMPORTANT:** measured without seat cushion!

## Quickie Easy 200

● = Standard    ○ = No cost option    □ = Cost option

Price in €

### SEAT UPHOLSTERY AND CUSHIONS

#### Seat Slings:

ESY020001	●	Seat sling with one utility bag	
ESY020004	○	Without seat sling	

#### Seat cushion black cover:

Foam composite: (medium soft)	Thickness:	ESY020008	□	30 mm	ESY020009	□	50 mm	
-------------------------------	------------	-----------	---	-------	-----------	---	-------	--

#### JAY Cushions - for complete product range see Jay order forms:

JAY020002	□	Jay Basic	
JAY020003	□	Jay Soft Combi P	
JAY020001	□	Jay Easy Visco Curved Cushion	

#### Seat accessories

ESY020016	□	Positioning belt with lock	
-----------	---	----------------------------	--

### BACKREST

**Back Height Range:** Choose back height factory setting (mm).

275 - 300	300 - 350	350 - 400	400 - 450	450 - 475
○ 275	○ 300	○ 350	○ 400	○ 450
○ 300	○ 325	○ 375	○ 425	○ 475
○ 350	○ 350	○ 400	○ 450	

Height adjustments are possible in the range: 275-300 mm; 300-350 mm; 350-400 mm, 400-450 mm and 450-475 mm.

ESY030031	○	Std back 90° (non angle adjustable)	
ESY030112	□	Angle Adjustable Back - choose factory setting (1)	(only with desk armrest & alu sideguard)
ESY030524	□	80°	
ESY030525	□	85°	
ESY030053	□	90°	
ESY030054	□	95°	
ESY030055	□	100°	
ESY030056	□	105°	(not in combination with desk sideguard)
ESY030011	□	Reclining back; -3° to +20° with gas spring	(Only with desk armrest - BH from 30 cm; not with F16)
		Lumbar bend	
ESY030012	□	5° frontwards	(from BH 350)
ESY030013	□	5° backwards	(from BH 350)

(1) Back angle 100°-105° not in combination with COG 25 mm and one arm drive; not with F16  
not with scalamobil frame

For more information  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
© Sunrise Medical 2025

## SEAT CUSHIONS



Seat cushion black cover



Seat cushion with water  
resistant / breathable cover



Medicin coat



Jay Basic  
(picture shown without upholstery)



Jay Soft Combi P  
(picture shown without upholstery)



Jay Easy Visco  
(picture shown without upholstery)

## SEAT ACCESORIES



Positioning belt with lock

## Quickie Easy 200

● = Standard

○ = No cost option

□ = Cost option

Price in €

### Pushhandles: Choose only one option

ESY030014	○	Integral Pushhandles long	
ESY030015	○	Integral Pushhandles short (until June 30th 2020)	(not in combination with drumbrake; R20)
ESY030016	○	Omit Pushhandles	(Only from back height 35 cm, not with drum brake / stabilizer bar / height adj push handles/ R20)
ESY030017	□	Fold-Down Pushhandles	(Only from back height 35 cm, not with drum brake / stabilizer bar / height adj push handles)
ESY030018	□	Height adjustable pushhandles	
ESY030020	□	Autofolding stabilizer bar	(Not with sw 32; 34 and 52; 54 cm / height adj or fold down push handles)
ESY030021	□	JAY Hardware height adjustable pushhandles	(Not with headrest EAY030022 & EAY030023)

### Head rest

ESY030022	□	medium	(Only with height adj. push handles)
ESY030023	□	large	(Only with height adj. push handles)

### Backrest upholstery

ESY030316	○	EXO Backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag, black (EXO Evo version)	
		Colour, binding tape: ● Black ○ Orange ○ Silver ○ Red ○ Blue	
ESY030317	□	EXO PRO Backr. sling, padded, adj., breathable spacer fabric, accessory bag, black (EXO Evo version)	
		Colour, binding tape: ● Black ○ Orange ○ Silver ○ Red ○ Blue	
ESY030300	○	Without backrest upholstery	

### JAY BACKS plus ACCESSORIES - for complete product range see Jay order forms

#### Jay 3 Back - Shallow Contour (SC)

Mid Thoracic:	<input type="checkbox"/>	JAY030010	Medium - 340 mm	<input type="checkbox"/>	JAY030033	Tall - 380 mm			
Upper Thoracic:	<input type="checkbox"/>	JAY030012	Short - 420 mm	<input type="checkbox"/>	JAY030019	Medium - 460 mm	<input type="checkbox"/>	JAY030025	Tall - 500 mm
Shoulder Height:	<input type="checkbox"/>	JAY030015	Short - 530 mm	<input type="checkbox"/>	JAY030022	Medium - 570 mm	<input type="checkbox"/>	JAY030028	Tall - 610 mm

#### Jay 3 Back - Mid Contour (MC)

Mid Thoracic:	<input type="checkbox"/>	JAY030011	Medium - 340 mm	<input type="checkbox"/>	JAY030018	Tall - 380 mm			
Upper Thoracic:	<input type="checkbox"/>	JAY030013	Short - 420 mm	<input type="checkbox"/>	JAY030020	Medium - 460 mm	<input type="checkbox"/>	JAY030026	Tall - 500 mm
Shoulder Height:	<input type="checkbox"/>	JAY030016	Short - 530 mm	<input type="checkbox"/>	JAY030023	Medium - 570 mm	<input type="checkbox"/>	JAY030029	Tall - 610 mm

#### JAY J3 Back - Mid Deep Contour (MDC)

Mid Thoracic:	□	JAY030040	Short - 30 cm	□	JAY030044	Medium - 340 mm	□	JAY030041	Tall - 380 mm	
Upper Thoracic:	□	JAY030045	Short - 420 mm		JAY030042	Medium - 460 mm	□	JAY030046	Tall - 500 mm	
Shoulder Height:	□	JAY030043	Short - 53 cm	□	JAY030047	Medium - 570 mm	□	JAY030048	Tall - 610 mm	

#### JAY J3 Back - Mid Deep Contour (MDC) (Extra Deep, 180 mm Contour Depth)

Mid Thoracic:	□	JAY030049	Tall - 380 mm	
Upper Thoracic:	□	JAY030050	Medium - 460 mm	

#### Jay 3 Back - Deep Contour (DC)

Upper Thoracic:	□	JAY030014	Short - 420 mm	□	JAY030021	Medium - 460 mm	□	JAY030027	Tall - 500 mm	
Shoulder Height:	□	JAY030017	Short - 530 mm	□	JAY030024	Medium - 570 mm	□	JAY030030	Tall - 610 mm	

#### Jay 3 Back Headrests

WHT030415	□	Whitmyer Headrest, COBRA XTRA, flip back	(only with JAY J3 backs)
WHT030416	□	Whitmyer Headrest, AXYS, Standard	(only with JAY J3 backs)

### ARMRESTS

(choose only one option, not in combination with sideguards)

ESY040001	□	Single Post Height Adjustable - Adult	
ESY040002	□	Length 250 mm	
	□	Length 360 mm	
	□	Desk armrest with long and short pad	(only with 0° & 1° camber)
		Height-adjustable	
ESY040003	□	Short 250 mm	
ESY040004	□	Long 310 mm	
		Fix height	
ESY040005	□	Short 250 mm	
ESY040006	□	Long 310 mm	
ESY040007	□	Fix height, without pad	
ESY040008	○	Quad release for height adjustable armrest	(Not possible with single post height)

#### Armrest accessories

ESY040016	□	Therapy tray SW 380 - 480 mm	(Only with armpad)
ESY040017	□	Therapy tray flip up SW 340 - 500 mm (assembled to the right side)	(Only with armpad)
ESY040018	□	Therapy tray flip up SW 340 - 500 mm (assembled to the left side)	(Only with armpad)

For more information

[www.SunriseMedical.eu](http://www.SunriseMedical.eu)

© Sunrise Medical 2025

## PUSHHANDLES



Integral push handles  
long



Fold down push  
handles



Height adjustable  
push handles for JAY

## JAY BACKS



*EXO backrest sling (EXO  
Evo version)*



*EXO PRO backrest sling  
(EXO Evo version)*

JAY J3 Back  
Shallow Contour



50 mm



JAY J3 Back  
Mid Contour



80 mm



JAY J3 Back  
Mid Deep Contour



150/180 mm



JAY J3 Back  
Deep Contour



150 mm



## ARMREST



Single Post Height Adj. -  
Adult



Desk Armrest

Quickie Easy 200			
<input checked="" type="radio"/> = Standard <input type="radio"/> = No cost option <input type="checkbox"/> = Cost option		Price in €	
<b>SIDEGUARDS</b>			
ESY040009	<input type="checkbox"/> Plastic regular, detachable	(Quickie style)	
<b>Aluminium sideguards</b> (Easy style)			
ESY040010	<input type="checkbox"/> without fender		
ESY040012	<input type="checkbox"/> with fender	(only in combination with alu sideguard easy style)	
ESY040015	<input type="checkbox"/> toolheight adjustable armrest	(only in combination with alu sideguard easy style, not with angle adj. back)	
ESY040013	<input type="checkbox"/> Aluminium sideguard streight without fender (1)	(Xenon style)	
ESY040014	<input type="checkbox"/> Aluminium sideguard streight with fender (1)	(Xenon style)	
ESY040104	<input type="checkbox"/> Carbon sideguards with fender (1)		
<b>FOOTREST &amp; FOOTPLATES</b>			
<b>Footplates</b>			
ESY050006	<input type="radio"/> Angle adjustable footplates (composite)	(pair; only from sw 360mm; LLL: 320-500mm)	
ESY050007	<input type="checkbox"/> Aluminium Footplate divided	(only from seat width 340mm; LLL: 220-510)	
ESY050008	<input type="checkbox"/> Aluminium Footplate platform	(not in combination with the ELR; only from seat width 340mm; LLL: 220-510)	
ESY050009	<input type="checkbox"/> Swing Away mechanism for footboard	(only with platform footrest alu; from LLL 390mm)	
ESY050062	<input type="checkbox"/> Locking system, platform footplate, aluminium		
<b>Platform footrest, angle adjustable</b> (only for seat width 320 to 500 mm)			
ESY050010	<input type="checkbox"/> Aluminium board (XENON performance)	(depth adjustable; LLL: 320-480mm)	
ESY050011	<input type="checkbox"/> Composite board	(depth adjustable; LLL: 320-480mm)	
ESY050033	<input type="checkbox"/> Carbon board	(depth adjustable; LLL: 320-480mm)	
ESY050012	<input type="checkbox"/> Foot positioning plate	(only in combination with Platform)	
<b>Autofolding</b> (only sw 320 to 500 mm)			
ESY050027	<input type="checkbox"/> Composite board	(angle & depth adjustable; LLL: high frame from 350 - 49= mm, low frame from 310 - 460 mm)	
ESY050028	<input type="checkbox"/> Carbon board	(angle & depth adjustable; LLL: high frame from 350 - 49= mm, low frame from 310 - 460 mm)	
<b>LOWER LEG LENGTH (LLL)</b> Attention: LLL has to be 250 mm lower than the front seat height			
<div> <input type="radio"/> 220    <input type="radio"/> 230    <input type="radio"/> 240    <input type="radio"/> 250    <input type="radio"/> 260    <input type="radio"/> 270    <input type="radio"/> 280    <input type="radio"/> 290  <input type="radio"/> 300    <input type="radio"/> 310    <input type="radio"/> 320    <input type="radio"/> 330    <input type="radio"/> 340    <input type="radio"/> 350    <input type="radio"/> 360    <input type="radio"/> 370  <input type="radio"/> 380    <input type="radio"/> 390    <input type="radio"/> 400    <input type="radio"/> 410    <input type="radio"/> 420    <input type="radio"/> 430    <input type="radio"/> 440    <input type="radio"/> 450  <input type="radio"/> 460    <input type="radio"/> 470    <input type="radio"/> 480    <input type="radio"/> 490    <input type="radio"/> 500    <input type="radio"/> 510 </div>			
<b>LLL between 220 mm to 350 mm only as high mounted version</b> All measurements in this table are in mm			
<b>FRAME</b>			
ESY010001	<input checked="" type="radio"/> Standard Fixed Frame (70°) <input type="radio"/> Short frame 540 mm (SD 340 - 420 mm) <input type="radio"/> Medium frame 580 mm (SD 340 - 460 mm) <input type="radio"/> Long frame 620 mm (SD 340 - 500 mm)	(Not with SD 440-500) (Not with SD 480-500)	
ESY010004	<input type="radio"/> Tie down labels (crash test)		
ESY070851	<input type="checkbox"/> Amputee axle position		
<b>e-mpulse R20</b>			
EMO000203	<input type="checkbox"/> e-mpulse R20 added	(Please use separate OF to order the R20.)	
EMO000204	<input type="checkbox"/> e-mpulse R20 mounted	(Please use separate OF to order the R20.)	
<b>EMPULSE F55 BUNDLING</b>			
EMO000101	<input type="radio"/> With F55 order, not mounted to the wheelchair	(Please use separate OF to order the F55.)	
EMO000102	<input type="radio"/> With F55 order, mounted to the wheelchair	(Please use separate OF to order the F55.)	
<b>EMPULSE F35 BUNDLING</b>			
EMO000301	<input type="radio"/> With F35 order, not mounted to the wheelchair	(Please use separate OF to order the F35.)	
EMO000302	<input type="radio"/> With F35 order, mounted to the wheelchair	(Please use separate OF to order the F35.)	
(1) <b>Effective seat width (measured between the sideguards) is seat width + 20 mm</b>		For more information <a href="http://www.SunriseMedical.eu">www.SunriseMedical.eu</a> © Sunrise Medical 2025	

## SIDEGUARDS



Plastic regular, detachable



Aluminium sideguard w/o fender



Aluminium sideguard with fender



Aluminium sideguard with fender



Carbon Seitenteil mit Kleiderschutz

## FOOTREST & FOOTPLATES



Angle adjustable footplates (composite)



Aluminium Footplate divided



Aluminium Footplate platform



Aluminium board (XENON performance)



Carbon board



Foot positioning plate



Autofolding carbon board

## EMPULSE F35 & F55 BUNDLING - NOGO MATRIX FOR 24" REAR WHEELS

COG	Seat Depth						
	380	400	420	440	460	480	500
25	✓	✓	✓	✓	✓	✓	✓
45	✓	✓	✓	✓	✓	✓	✓
65	X	✓	✓	✓	✓	✓	✓
85	X	X	✓	✓	✓	✓	✓

## EMPULSE F35 & F55 BUNDLING - NOGO MATRIX FOR 25" REAR WHEELS

COG	Seat Depth						
	380	400	420	440	460	480	500
25	✓	✓	✓	✓	✓	✓	✓
45	X	✓	✓	✓	✓	✓	✓
65	X	X	✓	✓	✓	✓	✓
85	X	X	X	✓	✓	✓	✓



Quickie Easy 200		● = Standard    ○ = No cost option    □ = Cost option	Price in €
<b>CENTRE OF GRAVITY (COG)</b>		(Measured from front of back post; choose only one option)	
ESY070250	<input type="radio"/> 25 mm		
ESY070450	<input type="radio"/> 45 mm		
ESY070650	<input type="radio"/> 65 mm		
ESY070850	<input type="radio"/> 85 mm		
<b>FRAME ACCESSORIES</b>		(Choose either tip assist or anti- tip tubes, not both)	
<b>Tip assist (step tube)</b>			
ESY090001	<input type="checkbox"/> Left		
ESY090002	<input type="checkbox"/> Right		
ESY090003	<input type="checkbox"/> Pair		
ESY090004	<input type="checkbox"/> Crutch holder with Velcro	(only in combination with tip assist)	
ESY090005	<input type="checkbox"/> Transit wheels		
<b>Anti-tip tubes swing-away</b>			
ESY090006	<input type="checkbox"/> Left		
ESY090007	<input type="checkbox"/> Right		
ESY090008	<input type="checkbox"/> Pair		
<b>PREPARATION FOR ATTITUDE ADD-ON BIKE</b>			
ESY090211	<input type="checkbox"/> Preparation for ATTITUDE, w/o ATTITUDE order (1)	(docking pins assembled) (please include Attitude OF)	
ESY090212	<input type="radio"/> Preparation for ATTITUDE, with ATTITUDE order (1)	(docking pins assembled; ATTITUDE order seperately) (please include Attitude OF)	
<b>FORKS</b>			
(compare with your choice on the first page)			
ESY080001	<input type="radio"/> Fork standard, black (Carbotecture for 111 mm)	(not with powder add ons)	
	<input type="checkbox"/> Orange <input type="checkbox"/> Silver <input type="checkbox"/> Red <input type="checkbox"/> Blue		
ESY080029	<input type="radio"/> Castor fork, Aluminium, black		
	ESY080042 <input type="checkbox"/> Colour Cast. Fork Silver		
ESY080002	<input type="checkbox"/> Fork with quick release axle for castors		
<b>CASTOR</b>			
(compare with your choice on the first page)			
ESY080003	<input type="radio"/> Solid	(for 4" & 5" only 45 mm fork version)	All caster sizes
ESY080004	<input type="radio"/> Soft solid	(for 45 mm fork version only)	(Only for 5" & 6"; not with Frog Leg)
ESY080008	<input type="radio"/> Pneumatic		(Only for 6" & 7" castors)
ESY080005	<input type="checkbox"/> Solid Light ups		(Not available for 6"/7")
<b>Castor with aluminium rim, silver</b>			
ESY080006	<input type="checkbox"/> Solid	(Only for 4" and 5")	103 €
ESY080007	<input type="checkbox"/> Soft roll	(Only for 4" and 5"; 45 mm fork only)	648 €
	Colour, rim: <input type="checkbox"/> Black <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue		53 €
<b>REAR WHEELS</b>			
(compare with your choice on the first page; choose only one wheel type)			
<b>NOTE: for users with higher user weight than 100 kg, we recommend cross spoked wheels</b>			
ESY070001	<input type="radio"/> Universal, cross spoked Wheels	(22"; 24"; 26")	
ESY070002	<input type="radio"/> Design Wheels	(24"; 26")	
ESY070003	<input type="radio"/> Lightweight Wheels	(22"/24"/25"; max. user weight: 125 kg)	
ESY070004	<input type="checkbox"/> Proton Wheels	(22"/24"/25"; max. user weight: 100 kg)	
ESY070005	<input type="checkbox"/> Spinergy Wheels	(24"; 25")	
ESY070006	<input type="checkbox"/> Rear wheels with drum brake	(22"; 24")	
ESY070007	<input type="checkbox"/> Mountain Bike wheels	(only with 24", mountain bike tyres & aluminium rim, black coated)	
ESY070008	<input type="checkbox"/> Additional mountain bike wheels pair	(only with 24"; comes always with mountainbike tyres )	
<b>One arm drive</b>		22" & 24"; only with stainless steel handrim and 0° camber	
ESY070009	<input type="checkbox"/> Left side		
ESY070010	<input type="checkbox"/> Right side		
ESY070011	<input type="checkbox"/> Omit Rear Wheels		
<b>TYRES</b>			
(choose only one tyre type)			
ESY070015	<input type="radio"/> Schwalbe Right Run	(Not in 25")	
ESY070107	<input type="checkbox"/> Schwalbe One	(Not in 22")	
ESY070018	<input type="radio"/> Full Poly Tire	Proton / Light Weight wheel only 24" & 25"; Spinergy wheels only with 24"; not with Mountainbike wheel)	
ESY070019	<input type="checkbox"/> Schwalbe Marathon Plus Evolution	(Not in 22")	
		For more information <a href="http://www.SunriseMedical.eu">www.SunriseMedical.eu</a> © Sunrise Medical 2025	

## FRAME ACCESSORIES



*Tip assist left*



*Anti-tip tubes SA*

## REAR WHEELS



*Universal Cross spoked*



*Design Wheel*



*Light weight wheel*



*Proton Wheel*



*Spinergy Wheel*



*Right Run*



*Schwalbe One*



*Marathon Plus  
(puncture proof)*



*Mountain bike rear  
wheel*

## Quickie Easy 200

● = Standard      ○ = No cost option      □ = Cost option

Price in €

### HANDRIM

ESY070030	○ Aluminium silver anodized	
ESY070031	○ Aluminium silver powdercoated	
ESY070032	□ Coloured aluminium pushrim. colour code: .....	
ESY070033	○ Stainless steel	
ESY070034	□ Titanium (not with 22"; 25" only in combination with Spinergy)	
ESY070035	□ Titanium clear powdercoated (not with 22"; 25" only in combination with Spinergy)	
ESY070036	□ Supergrip (22" only with universal and drum brake wheels)	
ESY070037	□ Maxgripp (only with 22" & 24", only handrim mounting wide)	
ESY070212	□ Handrims, Ellipse 3R, black (ergonomic handrims with rubber-strip; only 24" and 25")	
ESY070038	□ Handrim cover (not with 25", Supergrip, Maxgripp)	

### Handrim mounting

ESY070050	● Close (12 - 17mm) (not with Maxgripp)	
ESY070051	○ Wide (27mm) (Not with light weight, Proton & Spinergy wheel and alu handrim)	
ESY070052	□ Extra Close (7mm) (only with stainless steel, only with 24" Design & Universal wheel)	

not in combination with Spinergy, Proton and Lightweight Wheels and stainless steel Handrim

### WHEEL CAMBER

(choose only one option)

ESY070053	● 0°	
ESY070054	□ 1° (Not in combination with one arm drive)	
ESY070055	□ 2° (Not in combination with desk sideguard; one arm drive)	
ESY070056	□ 3° (Not in combination with desk sideguard; one arm drive)	
ESY070057	□ 4° (Not in combination with desk sideguard; one arm drive)	

### SPOKE GUARDS (only available in 22" & 24")

(Not with Spinergy Wheels)

ESY070060	□ Transparent	
	ESY070073 □ Galactica	
ESY070063	□ Princess	ESY070085 □ Jungle Fever
ESY070064	□ Happy Lions	ESY070087 □ Feather
ESY070065	□ Monkeys	ESY070076 □ Pirates
ESY070066	□ Oriental	ESY070089 □ Chrome wheel
ESY070067	□ Superhero	ESY070078 □ Graffiti
ESY070068	□ Hot wheels	ESY070079 □ Basketball
ESY070069	□ Pool	ESY070080 □ Flying water
ESY070070	□ Silk	ESY070091 □ Space
ESY070071	□ Tiger	ESY070082 □ Zebra

### WHEEL LOCKS (Choose only one option)

ESY060001	○ Standard wheel lock	
ESY060002	□ Knee lever wheel lock	
EAY060023	□ Compact wheel lock EVO	
ESY060004	□ Extension for knee lever wheel lock, slip-on, swing-away (Only on kneelever wheel lock)	
<b>One arm wheel lock</b>		
ESY060005	□ Left side	
ESY060006	□ Right side	

### ACCESSORIES

ESY090010	● Tool kit	
ESY090011	□ High pressure pump	
ESY090012	□ Air pump	
ESY090013	□ Back pack	
ESY090016	□ Mobile phone bag	

(1) Only with forks for power add on, not with SB 320, 480-54 0mm, not with SH front 430-490 mm, not with reclining an angle adj back, transit wheels, anti tips an tip assist)

For more information  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
 © Sunrise Medical 2025

## HANDRIMS



*Max Grepp*



*Supergrip*



*Ellipse 3R*

## WHEEL LOCK



*Standard Wheel lock*



*Compact wheel lock*



*Knee lever wheel lock*



*One arm wheel lock right*

## ACCESSORIES



*Tool kit*



*High pressure pump*




*Air pump*



*Back pack*



*Mobile phone pocket*

Quickie Easy 200							
<input checked="" type="radio"/> = Standard <input type="radio"/> = No cost option <input type="checkbox"/> = Cost option				Price in €			
<b>COLOURS <sup>(1)</sup></b> Please consult our colour catalogue to review your colour choice							
	Frame	Handrim	Hub/rim		Frame	Handrim	Hub/rim
02 Brilliant Silver	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	27 Graphitemetallic	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
03 Silver Matt	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	29 Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
06 Sun Yellow	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	31 Matt Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
14 Blue	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	32 Signal White	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
15 Brilliant Marine	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	33 Pink	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
16 Matt Marine	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	34 Dark Chocolate Metallic	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
19 Brilliant Purple	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	35 Skyblue	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
21 Red Ruby	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	36 Cyan-Metallic	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22 Brilliant Red	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	37 Marsred (Orange)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
23 Matt Red	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	40 Kiwi	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24 Tutti Frutti Orange	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	41 Flamingo-Glitter	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
				42 Matt-Anthracite	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
				46 Mamba Green (new)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
(1) Handrim: only with handrim colour powder coated; Hub/rim: only for universal wheel							
<b>STICKER COLOURS</b>							
ESY100550	<input type="radio"/>	Sticker colour - black					
ESY100551	<input type="radio"/>	Sticker colour - silver					
<b>Anodized colour axle plate</b>							
ESY100580	<input checked="" type="radio"/>	Black					
ESY100581	<input type="checkbox"/>	Silver					
ESY100582	<input type="checkbox"/>	orange					
<b>24 " LIGHT WEIGHT WHEEL COLOURS (anodized)</b>							
	<input checked="" type="radio"/>	Black					
ESY100336	<input type="checkbox"/>	Anodized, orange					
ESY100337	<input type="checkbox"/>	Anodized, blue					
ESY100338	<input type="checkbox"/>	Anodized, red					
<b>Comments for special build:</b>							
<b>Examples of available build 4 me options</b> Abduction frame fixed (Build 4 ME # 0762) Back height 500-550 mm (Build 4 ME # 0320) Front SH 550 mm (Build 4 ME # 0823) 4"/5" solid castors with alu rim red or blue anodized (FMG-MG01-2994-01)							
<div style="text-align: right;">  </div>							
<b>Easy 200 - Base Price</b>							
Option price							
Packaging and transport							
Exclusive VAT							
VAT.							
<b>Total Price</b>							
Sunrise Medical Customer service Address: Ul. Walczaka 114a,66 400, Gorsaw, Gorzów Wlkp. Polen Tel. +48 95 737 25 14							
E-mail: <a href="mailto:cs.pl@sunrise-medical.pl">cs.pl@sunrise-medical.pl</a> Website: <a href="http://www.sunrisemedical.eu">www.sunrisemedical.eu</a>							
<a href="http://www.SunriseMedical.eu">www.SunriseMedical.eu</a> © Sunrise Medical 2025							

## ANODIZED COLOUR AXLE PLATE



Axle plate  
anodized, orange



Axle plate  
anodized, silver



Axle plate  
anodized, black

## 24 " LIGHT WEIGHT WHEEL COLOURS



*Lightweight wheel, anodized orange*



*Lightweight wheel, anodized blue*



*Lightweight wheel, anodized red*