



WHERE YOU WANT TO BE.

# LUXURY IN MOTION

**S-SERIES**  
— PREMIUM SCOOTER RANGE —

**S400**

**S425**

**S700**



# INNOVATIVE S-SERIES SCOOTERS

---

From the shopping mall to the country, the new luxury S-Series by STERLING are ready for any journey. Packed full of brilliant features with superb driving dynamics, and looking sharp from every angle, there's an S-Series mobility scooter to suit everyone.

*S425*










# ADVANCED LED LIGHTING

---

**Powerful all-around LED lighting, approved for road use. With lifetime performance, changing bulbs is a thing of the past.**





# INNOVATION IN DESIGN

---

With advanced LED lighting technology, the S-Series scooters are 400 times more energy efficient than standard bulbs, meaning more battery power and range every trip.

---

**Safety as standard, S-Series models all feature rear high visibility LED indicators, position and brake lights.**



→ **Innovative active anti-tip wheels**  
Move up as you ride over obstacles,  
and keep you safe on slopes.





# S400

## COMPACT SCOOTER

---

The S400 **6kph** luxury compact scooter with big features, is highly maneuverable around shopping centres or busy high streets, and with the all-around suspension as standard makes every journey a comfortable one.

With **38 Ah** batteries the S400 is ideal for trips around town, plus with the bright, energy efficient LED's lighting your way, you're just as visible on those late-night shopping trips as you are during the day.





Max Speed: **6 kph**



Colour: **Silver Metallic**



Batteries: **38 Ah**



Wheel size: **25 cm**



Max. user weight: **136 kg**





Waterproof soft touch controls featuring easy to use speed adjustment, LED battery gauge and a powerful horn.


Unique wig-wag design means natural driving using finger or thumb control. The large ergonomically designed delta handlebar makes every journey a comfortable one.



# CONTROLS DESIGNED FOR YOU

---

Finding your perfect driving position is a breeze with the infinitely adjustable tiller. Clever design means the headlight automatically tilts as the tiller is adjusted.



Safety and style as standard, every S-Series scooter comes with two carbon effect mirrors

The removable basket with carry handle, useful for those shopping trips.

# S425

## COMPACT OUTDOOR SCOOTER

---

The luxury S425 crams in all of the features of its little brother, but with added power!

Taking the bright, front/rear LED lights, comfortable all-around suspension and combining it with a maximum speed of **12 kph**, the S425 is suitable for road use as well as on pavements.

The increased wheel height will help to tackle those usual outdoor obstacles, and longer journeys can be handled with ease as the larger **55 Ah** batteries will help to power your way.







Max Speed: **12 kph**



Colour: **Carbon Metallic**



Batteries: **55 Ah**



Wheel size: **30 cm**



Max. user weight: **150 kg**



# LUXURY SEATING

---

**The multi-adjustable luxury seating is designed to make every journey a comfortable one.**





# COMFORT ON EVERY TRIP

---

The S-SERIES advanced seating solution, takes comfort to the highest level, from the seat height, depth and recline adjustments, to the flip up, width, angle and depth adjustable comfort armrests – every feature has been designed for comfort – for you.



→ Infinite seat depth adjustable



→ Quick backrest angle adjustment



→ Seat rotation enables comfortable transfers

# S700

## PERFORMANCE OUTDOOR SCOOTERS

---

The luxury S700, with superior outdoor performance is your perfect companion for outdoor adventures, to take you where you want to go.

There's no need to worry about distance either as the large **75 Ah** battery size will power an extended range and a lively top speed of **15 kph**.

The high performance design means gradients, and higher kerbs are tackled with ease and the large tyres and advanced 4-wheel suspension provides a smooth, rewarding drive each and every time.







Max Speed: **6 kph / 12 kph / 15 kph**



Colour: **Silver- / Carbon- / Black Metallic**



Batteries: **75 Ah**



Wheel size: **33 cm**



Max. user weight: **160 kg**

# TECHNICAL DATA

S-SERIES	S400	S425	S700
Max speed	6 kph	12 kph	6 kph / 12 kph / 15 kph
Max. user weight	136 kg	150 kg	160 kg
Battery size	38 Ah x 2	55 Ah x 2	75 Ah x 2
Charger	5 A	8 A	8 A
Max. potential range*	36 km	42 km	55 km
Max. power (peak)	1000 W	1200 W	2700 W
All around LED lighting	Yes	Yes	Yes
Easy adjust tiller & wig wag control	Yes	Yes	Yes
Suspension	Yes	Yes	Yes - advanced
Tyres	25 cm pneumatic	30 cm pneumatic	33 cm pneumatic
Overall width / length	64 cm / 134 cm	66 cm / 135 cm	67 cm / 150 cm
Seat width x depth	46 cm x 46 cm	49 cm x 49 cm	49 cm x 49 cm
Seat height range	45 cm - 50 cm	45 cm - 50 cm	45 cm - 50 cm
Minimum turning distance	1.9 m	2.0 m	2.2 m
Colour	Silver Metallic	Carbon Metallic	Silver / Black / Carbon Metallic
Scooter weight inc. batteries	109 kg	117 kg	146 kg

\* Max potential range based on industry benchmark range figures for similar products and battery sizes

All over Europe, thousands of people have made their lives easier with a STERLING mobility scooter by Sunrise Medical. For over 30 years, STERLING have designed and built mobility scooters to the highest standards, giving the comfort and safety every customer needs. For further information on product options please check the order form on our website. You can view this document also in a larger font on a PDF-file at [www.SunriseMedical.eu](http://www.SunriseMedical.eu)



## Sunrise Medical GmbH

Kahlbachring 2-4  
 69254 Malsch/HD  
 Tel.: +49 (0) 72 53 / 9 80-0  
 Fax: +49 (0) 72 53 / 9 80-222  
 E-Mail: [international@sunrisemedical.de](mailto:international@sunrisemedical.de)  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)

